

Hanover Prosthetic Fabrication  
 219 Hillsborough Rd. Unit A  
 Carrboro, NC 27510  
 (336)260-2208  
 hanoverprofab.com  
 contact@hanoverprofab.com



## AK Test Socket Order Form

Hanover Acct# \_\_\_\_\_ Date : \_\_\_\_\_

Company: \_\_\_\_\_ Ship to: \_\_\_\_\_


\_\_\_\_\_  
 \_\_\_\_\_

Ordering Practitioner: \_\_\_\_\_  
 Email: \_\_\_\_\_ Contact Phone : \_\_\_\_\_


Patient last name or other I.D. : \_\_\_\_\_

PO # \_\_\_\_\_

**K-Level :** \_\_\_\_\_ **Amp. Level:** AK KD **Side:** Left Right **Patient Weight :** \_\_\_\_\_



mm	Loose	Tight
1		
2		
3		
4		
5		
6		
7		
IT to Distal _____ mm		




Liner Fitted \_\_\_\_\_  
 \_\_\_\_\_

Additional reduction over measurements \_\_\_\_\_ Ply  
 \_\_\_\_\_ %

**Test Socket Material:**  
 PETG-Std.  
 Other \_\_\_\_\_  
 \_\_\_\_\_

**Mounting Plate Alignment**



**Options:**

- Add distal mounting plate (4 hole Std.) - Other \_\_\_\_\_
- Standard AK expulsion valve - Other \_\_\_\_\_
- Lanyard base (KISS Delrin Lanyard Base Std.) - Other \_\_\_\_\_
- Shuttle lock (Fillauer 125232MM Std.) - Other \_\_\_\_\_

**Notes:** \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

3-5 deg. ADuct - Std.    5-7 deg. Flexion - Std.  
 Other \_\_\_\_\_ deg.    Other \_\_\_\_\_ deg.

**Knee unit used :** \_\_\_\_\_  
 \_\_\_\_\_